

# ChromTech Double Beam UV-6 Series Spectrophotometer

**UV-6 Series** are advanced double beam design consisting of four models:

Stand-alone models: **UV6100** with 1.8nm fixed bandwidth and **UV6300** with 1.0nm fixed bandwidth;

PC models: **UV-6100PC** with 1.8nm fixed bandwidth and **UV-6300PC** with 1.0nm fixed bandwidth.

Other specifications of the four styles are almost the same except bandwidth.

The two detectors are measuring sample and reference respectively and simultaneously for optimizing measurement accuracy.

They provide excellent performance for measurements in the range of 190nm

to 1100nm. They are suitable for pharmaceutical, biochemical and clinical lab

applications as well as routine applications such as quantitative analyses, kinetics, spectrum scanning, multiple components and DNA/Protein, PC Windows application software make these instruments versatile. All instruments provide excellent performance for measurements.



- To Stand-alone models, All software methods are included as built-in standard, thus eliminating the need for software options.
- Online software upgrade via internet helps to keep your software up-to-date.
- Data Download-to-PC software expands the data storage to unlimited.

**Stand-alone models of UV-6 Series have the same functions as UV-3 series have, see next page for details.**

Model	UV-6100	UV-6100PC	UV-6300	UV-6300PC
Wavelength Range	190-1100nm		190-1100nm	
Spectral Bandwidth	1.8nm	1.8nm	1.0nm	1.0nm
Optical System	Double Beam, Grating 1200 lines/mm			
Wavelength Accuracy	±0.3nm			
Wavelength Repeatability	0.2nm			
Scanning Speed	Hi, Med., Low. Max.3000nm/min			
Photometric Accuracy	±0.2% T			
Photometric Repeatability	±0.15% T			
Photometric Range	-0.3-3A, 0-200% T, 0-9999C conc.			
Stray Light	0.05% T			
Stability	±0.0005A/h			
Display	LCD (320X240)	PC Model	LCD (320x240)	PC Model
Baseline Flatness	±0.0005A (200-1000nm)			
Sample Compartment	Accommodates 100mm pathlength cuvette with optional holder			
Light Source	Halogen & Deuterium lamp (pre-aligned)			
Output	USB Port & Parallel Port (Printer)			
Power Requirement	AC 220V/50Hz or AC 110V/60Hz			
Dimensions (W x D x H)	600 x 450 x 200mm			
Weight	22kg		22kg	

The PC models come standard with Windows® based application software.


**PC models of UV-6 Series have the same functions as UV-3 series have, see next page for details.**




# ChromTech Local-Control Software For UV-6 Series

All software methods are included as built-in standard, thus eliminating the need for software options. Online software upgrade via Internet helps to keep your software up-to-date.

The local control software include functions as: Basic Mode, Quantitative, Wavelength Scan, Kinetics, DNA/Protein, Multi-wavelength Test and System Utilities.




## Main Menu




## Standard Curve

Up to 10 standard solutions may be used to establish calibration equation curve. There is a choice of four methods for fitting a curve through the calibration points: Linear fit, Linear fit through zero, square fit and cubic fit.



## Wavelength Scan

The wavelength scan intervals are 0.1, 0.2, 0.5, 1, 2, 5nm, and Hi, Medium and Low scan speeds are available. Scan speeds vary from 100 to 1000 nm/min. Wavelengths are scanned from high to low so that the instrument waits at high wavelength. This minimizes the degradation of UV sensitive samples. Precise control of filter and lamp changes means that their effects are not seen on the final scan. Post-run manipulation includes re-scaling axes, curve tracking and peak picking



## Kinetics

This mode may be used for time course scanning or reaction rate calculations. Abs. vs. time graphs are displayed on the screen in real time. Wait time and measurement time up to 12 hours may be entered with time intervals of 0.5, 1, 2, 5, 10, 30 seconds and 1 min. Post-run manipulation includes re-scaling, curve tracking and selection of the part of the curve required for the rate calculation. Rate is calculated using a linear regression algorithm before multiplying by the entered factor.

No	WL (nm)	Abs
1	500.0	0.87
	400.0	0.42
	300.0	0.81

## Multi-Wavelength

Up to 10 wavelengths may be entered, allowing the measurement of multiple wavelengths on a series of Samples.

No	Items	Result	Unit
1	A1	2.947	Abs
	A2	2.842	Abs
	Aref	0.638	Abs
	C-DNA	65.91	mg/mL
	C-Pro	1672	mg/mL
	Ratio	1.048	

## DNA/Protein Test

Concentration and DNA purity are calculated:

Absorbance ratios 260nm/280nm or 260nm/230nm

With optional subtracted absorbance at 320nm

DNA Concentration=62.9 x A260 - 36.0 x A280

or 49.1x A260 - 3.48x A230

Protein Concentration=1552 x A260 - 757.3 x A280

or 183 x A260 - 75.8 x A230


Other wavelengths and ratio factors may be entered.

# ChromTech PC- Control Software For UV-6 Series

The **ChromTech** Windows-based PC application software takes the best features of the stand-alone version plus more powerful data processing and expanded data collecting and storage capability. It comes standard with **ChromTech** PC models and is optional to stand-alone models.

## The PC application software offers:

1. Basic Photometric Mode
2. Quantitative test (standard curve)
3. Wavelength Scanning
4. Kinetics
5. DNA/Protein
6. Multi-wavelength Test
7. System Utility



## Kinetics(Abs vs. Time)

The Kinetics mode may be used for time course scanning or reaction rate calculations. Abs. vs. time graphs is displayed on the screen in real time. Wait time, measurement time and time intervals may be entered. Post-run manipulation includes re-scaling, curve tracking and selection of the part of the curve required for the rate calculation. Rate is calculated using a linear regression algorithm before multiplying by the entered factor.

Sample Name	260nm	280nm	C-DNA	Protein Ratio
Sample-1	4.1004	4.5957	0.4760179	7289454.3120
Sample-2	0.1390	0.2176	0.5559	4.2943 39.1040
Sample-3	0.1587	0.3123	0.5045	4.2199 39.2333
Sample-4	4.4217	4.5885	0.4710179	6475454.0293
Sample-5				
Sample-6				
Sample-7				
Sample-8				
Sample-9				
Sample-10				
Sample-11				
Sample-12				
Sample-13				
Sample-14				
Sample-15				
Sample-16				
Sample-17				
Sample-18				
Sample-19				
Sample-20				
Sample-21				
Sample-22				
Sample-23				
Sample-24				

## DNA/Protein


Concentration and DNA purity are quickly and easily calculated:

Absorbance ratios 260nm/280nm with optional subtracted absorbance at 320nm.

DNA concentration= 62.9 x A260 - 36.0 x A280


Protein concentration=1552 x A260 - 757.3 x A280

Other wavelengths and factors may be entered.



## Multi-wavelength


Up to 32 wavelengths can be selected and multiple samples can be measured.(Auto cell changer is required to run multiple samples automatically)



## Quantitative Test (Standard Curve)

Use up to 32 standards to establish standard curve. Four methods for fitting a curve:

1. Linear fit
2. Linear through zero
3. Square fit
4. Cubic fit



## Wavelength Scanning

Automatically record peaks and valleys. The quantity of channels is unlimited, you can simultaneously store curves as many as you want.

Post-run manipulation and processing includes:

1. Re-scaling axes, curve
2. 1<sup>st</sup> to 4<sup>th</sup> derivative
3. Smoothing, combination, zooming, overlap...