

毛細電泳層析管 CE Capillary Tubing



FLEXIBLE FUSED SILICA CAPILLARY TUBING

目錄編號	產品名稱	管柱內徑(μ)	管柱外徑(μ)
2000001	TSP002150	002	150
2000002	TSP005150	005	150
2000003	TSP005375	005	363
2000004	TSP010150	010	150
2000005	TSP010375	010	363
2000006	TSP015150	015	150
2000007	TSP015375	015	363
2000381	TSP020090	020	090
2000008	TSP020150	020	150
2000009	TSP020375	020	363
2001145	TSP025150	025	150
2000011	TSP025375	025	363
2000012	TSP030150	030	150
2000013	TSP030375	030	363
2000596	TSP040105	040	105
2000014	TSP040150	040	150
2000383	TSP040375	040	363
2000015	TSP050150	050	150
2000016	TSP050192	050	186
2000017	TSP050375	050	363
2000018	TSP075150	075	150
2000133	TSP075200	075	193
2000019	TSP075375	075	363
2000020	TSP100170	100	164
2000021	TSP100200	100	193
2000022	TSP100245	100	238
2000023	TSP100375	100	363
2000024	TSP150375	150	363
2000025	TSP180350	180	340
2000204	TSP200350	200	360
2000026	TSP250350	250	360
2000027	TSP320450	320	435
2000000	TSP530660	542	673
2000028	TSP530700	530	700
2000029	TSP700850	700	850

* 操作最高耐溫至 350°C

* 最低購買數量 10 公尺