


標準液槽 Standard Cells

標準石英液槽 (兩面磨光,兩面磨砂)

※ 適用波長：190nm-2500nm

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 標準石英液槽(含蓋)	45 x 12.5 x 3.5	1	10	0.35	Q11001	2 支/盒
	45 x 12.5 x 4.5	2	10	0.70	Q11002	2 支/盒
	45 x 12.5 x 7.5	5	10	1.7	Q11005	2 支/盒
	45 x 12.5 x 12.5	10	10	3.5	Q11010	2 支/盒
	45 x 12.5 x 22.5	20	10	7.0	Q11020	2 支/盒
	45 x 12.5 x 32.5	30	10	10.5	Q11030	2 支/盒
	45 x 12.5 x 42.5	40	10	14	Q11040	2 支/盒
	45 x 12.5 x 52.5	50	10	17.5	Q11050	2 支/盒
 標準石英液槽(含塞)	45 x 12.5 x 7.5	5	10	1.7	Q12005	2 支/盒
	45 x 12.5 x 12.5	10	10	3.5	Q12010	2 支/盒
	45 x 12.5 x 22.5	20	10	7.0	Q12020	2 支/盒
	45 x 12.5 x 32.5	30	10	10.5	Q12030	2 支/盒
	45 x 12.5 x 42.5	40	10	14.0	Q12040	2 支/盒
	45 x 12.5 x 52.5	50	10	17.5	Q12050	2 支/盒
	45 x 12.5 x 102.5	100	10	35.0	Q11100	2 支/盒
	45 x 12.5 x 102.5	100	10	35.0	Q12100	2 支/盒

標準玻璃液槽 (兩面磨光,兩面磨砂)

※ 適用波長：340nm-2500nm



產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 標準玻璃液槽(含蓋)	45 x 12.5 x 3.5	1	10	0.35	G11001	4 支/盒
	45 x 12.5 x 4.5	2	10	0.70	G11002	4 支/盒
	45 x 12.5 x 7.5	5	10	1.7	G11005	4 支/盒
	45 x 12.5 x 12.5	10	10	3.5	G11010	4 支/盒
	45 x 12.5 x 22.5	20	10	7.0	G11020	4 支/盒
	45 x 12.5 x 32.5	30	10	10.5	G11030	4 支/盒
	45 x 12.5 x 42.5	40	10	14	G11040	4 支/盒
	45 x 12.5 x 52.5	50	10	17.5	G11050	4 支/盒
 標準玻璃液槽(含塞)	45 x 12.5 x 7.5	5	10	1.7	G12005	4 支/盒
	45 x 12.5 x 12.5	10	10	3.5	G12010	4 支/盒
	45 x 12.5 x 22.5	20	10	7.0	G12020	4 支/盒
	45 x 12.5 x 32.5	30	10	10.5	G12030	4 支/盒
	45 x 12.5 x 42.5	40	10	14.0	G12040	4 支/盒
	45 x 12.5 x 52.5	50	10	17.5	G12050	4 支/盒
	45 x 12.5 x 102.5	100	10	35.0	G11100	4 支/盒
	45 x 12.5 x 102.5	100	10	35.0	G12100	4 支/盒

螢光液槽 Fluorometer Cells

螢光石英液槽 (四面磨光)



產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 螢光石英液槽(含蓋)	45 x 12.5 x 7.5	5	10	1.7	Q21005	2 支/盒
	45 x 12.5 x 12.5	10	10	3.5	Q21010	2 支/盒
	45 x 12.5 x 22.5	20	10	7.0	Q21020	2 支/盒
	45 x 12.5 x 32.5	30	10	10.5	Q21030	2 支/盒
	45 x 12.5 x 42.5	40	10	14	Q21040	2 支/盒
	45 x 12.5 x 52.5	50	10	17.5	Q21050	2 支/盒
	45 x 12.5 x 7.5	5	10	1.7	Q21005	2 支/盒
 螢光石英液槽(含塞)	45 x 12.5 x 7.5	5	10	1.7	Q22005	2 支/盒
	45 x 12.5 x 12.5	10	10	3.5	Q22010	2 支/盒
	45 x 12.5 x 22.5	20	10	7.0	Q22020	2 支/盒
	45 x 12.5 x 32.5	30	10	10.5	Q22030	2 支/盒
	45 x 12.5 x 42.5	40	10	14.0	Q22040	2 支/盒
	45 x 12.5 x 52.5	50	10	17.5	Q22050	2 支/盒
	45 x 12.5 x 7.5	5	10	1.7	Q22005	2 支/盒

螢光微量石英液槽 (四面磨光)


產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 螢光微量石英液槽(含蓋)	45 x 12.5 x 12.5	10	2	0.7	Q23102	2 支/盒
	45 x 12.5 x 12.5	10	4	1.4	Q23104	2 支/盒
 螢光微量石英液槽(含塞)	45 x 12.5 x 12.5	10	2	0.7	Q24102	2 支/盒
	45 x 12.5 x 12.5	10	4	1.4	Q24104	2 支/盒

微量液槽 Micro Cells

微量石英液槽 (黑邊避光)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 黑邊微量石英液槽(含蓋)	45 x 12.5 x 12.5	10	1	0.35	Q31101	2 支/盒
	45 x 12.5 x 12.5	10	2	0.7	Q31102	2 支/盒
	45 x 12.5 x 12.5	10	3	1	Q31103	2 支/盒
	45 x 12.5 x 12.5	10	4	1.4	Q31104	2 支/盒
	45 x 12.5 x 12.5	10	5	1.75	Q31105	2 支/盒
 黑邊微量石英液槽(含塞)	45 x 12.5 x 12.5	10	2	0.7	Q32102	2 支/盒
	45 x 12.5 x 12.5	10	4	1.4	Q32104	2 支/盒

微量石英液槽 (磨砂白邊)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 白邊微量石英液槽(含蓋)	45 x 12.5 x 12.5	10	1	0.35	Q41101	2 支/盒
	45 x 12.5 x 12.5	10	2	0.70	Q41102	2 支/盒
	45 x 12.5 x 12.5	10	3	1	Q41103	2 支/盒
	45 x 12.5 x 12.5	10	4	1.4	Q41104	2 支/盒
	45 x 12.5 x 12.5	10	5	1.7	Q41105	2 支/盒
 白邊微量石英液槽(含塞)	45 x 12.5 x 12.5	10	2	0.7	Q42102	2 支/盒
	45 x 12.5 x 12.5	10	4	1.4	Q42104	2 支/盒

避光流動液槽 Self Masking Flow Cells

流動石英液槽 (黑邊避光)

產品圖形及說明	外部尺寸(HxWxD) mm	光徑長 mm	光徑寬 mm	容積 ml	產品編號	包裝
 避光流動石英液槽	45 x 12.5 x 12.5	10	φ1.5	0.018	Q51105	2 支/盒
	45 x 12.5 x 12.5	10	φ2	0.03	Q51102	2 支/盒
	45 x 12.5 x 12.5	10	φ3	0.07	Q51103	2 支/盒
 避光流動石英液槽	45 x 12.5 x 7.5	5	4x12	0.24	Q52054	2 支/盒
	45 x 12.5 x 12.5	10	4x12	0.48	Q52104	2 支/盒
	45 x 12.5 x 22.5	20	4x12	1.0	Q52204	2 支/盒
	45 x 12.5 x 32.5	30	4x12	1.5	Q52304	2 支/盒
 避光流動石英液槽	45 x 12.5 x 12.5	10	φ1	0.01	Q53101	2 支/盒
	45 x 12.5 x 12.5	10	φ1.5	0.018	Q53105	2 支/盒
	45 x 12.5 x 12.5	10	φ2	0.03	Q53102	2 支/盒

其他液槽



Micro cell with black walls



Cylindrical cell



Large cell



Micro cell with frosted walls