



# 可調式微量吸管 Adjustable Volume Pipettors

## ChromTech® 單道可調式微量吸管

### Single-Channel Adjustable Volume Pipettors

芬蘭 Biohit OEM

| 產品圖形及說明  | 編號     | 產品型號           | 容量               | 單位增量         | 適用管嘴         |
|--|--------|----------------|------------------|--------------|--------------|
|  <p>單道可調式微量吸管</p> | 720005 | CTYD0.1~2.5    | 0.1-2.5 $\mu$ l  | 0.05 $\mu$ l | 10 $\mu$ l   |
|  | 720000 | CTYD0.5~10     | 0.5-10 $\mu$ l   | 0.1 $\mu$ l  | 10 $\mu$ l   |
|  | 720080 | CTYD2~20       | 2-20 $\mu$ l     | 0.5 $\mu$ l  | 300 $\mu$ l  |
|  | 720050 | CTYD10~100     | 10-100 $\mu$ l   | 1 $\mu$ l    | 300 $\mu$ l  |
|  | 720070 | CTYD20~200     | 20-200 $\mu$ l   | 1 $\mu$ l    | 300 $\mu$ l  |
|  | 720060 | CTYD100~1000   | 100-1000 $\mu$ l | 5 $\mu$ l    | 1000 $\mu$ l |
|  | 720110 | CTYD1000~5000  | 1-5 ml           | 50 $\mu$ l   | 5 ml         |
|  <p>可調式微量吸管架</p> | 721001 | Pipette Holder | 5 支              |              |              |
|  |        |                |                  |              |              |
|  |        |                |                  |              |              |
|  |        |                |                  |              |              |
|  |        |                |                  |              |              |

## 多道可調式微量吸管

### Multi-Channel Adjustable Volume Pipettors

芬蘭 Biohit OEM

| 產品圖形及說明   | 編號     | 產品型號   | 容量                  | 單位增量        | 適用管嘴        |
|---|--------|--------|---------------------|-------------|-------------|
|  <p>八道可調式微量吸管</p>  | 720210 | CTYM8  | 8 道 0.5-10 $\mu$ l  | 0.1 $\mu$ l | 10 $\mu$ l  |
|   | 720220 | CTYM8  | 8 道 5-50 $\mu$ l    | 0.5 $\mu$ l | 300 $\mu$ l |
|   | 720240 | CTYM8  | 8 道 50-300 $\mu$ l  | 5 $\mu$ l   | 300 $\mu$ l |
|   |        |        |                     |             |             |
|  <p>十二道可調式微量吸管</p> | 720310 | CTYM12 | 12 道 0.5-10 $\mu$ l | 0.1 $\mu$ l | 10 $\mu$ l  |
|   | 720320 | CTYM12 | 12 道 5-50 $\mu$ l   | 0.5 $\mu$ l | 300 $\mu$ l |
|   | 720340 | CTYM12 | 12 道 50-300 $\mu$ l | 5 $\mu$ l   | 300 $\mu$ l |
|   |        |        |                     |             |             |