


固相萃取管 SPE Tube



Chromabond®

Phase	Column Material	1 ml		3 ml		6 ml	
		50 mg(100/PK)	100 mg(100/PK)	200 mg(50/PK)	500 mg(50/PK)	500 mg(30/PK)	1000 mg(30/PK)
C18	PP		730001	730002	730003	730004	730005
	Glass		730001G	730002G	730003G	730004G	730005G
C18 f (fast flow)	PP				730008		730009
	Glass				730008G		730009G
C18 ec (endcapped)	PP		730011	730012	730013	730014	730015
	Glass		730011G	730012G	730013G	730014G	730015G
C18 ec f (endcapped, fast)	PP				730018		730010
	Glass				730018G		730010G
C18 Hydra	PP	730294	730295	730296	730297	730299	730300
	Glass	730294G	730295G	730296G	730297G	730299G	730300G
C8	PP		730021	730022	730023	730024	
	Glass		730021G	730022G	730023G	730024G	
C4	PP		730225	730226	730227		
	Glass		730225G	730226G	730227G		
C2	PP		730169	730220	730221		
	Glass		730169G	730220G	730221G		
C ₆ H ₁₁ ec	PP		730440	730441	730442	730443	730444
	Glass		730440G	730441G	730442G	730443G	730444G
C ₆ H ₅	PP		730083		730084		
	Glass		730083G		730084G		

固相萃取管 SPE Tube

Chromabond®

NH₂ ↔Aminopropyl **PCA** ↔Propylcarboxylic acid **Alox N** ↔Aluminium oxide
DMA ↔Dimethylamino **PSA** ↔Propylsulphonic acid **PA** ↔Polyamide
CN ↔Cyanopropyl **SA** ↔Cation exchanger **HR-P** ↔Polystyrene
NO₂ ↔Nitrophenyl **SB** ↔Anion exchanger

Phase	Column Material	1 ml		3 ml		6 ml	
		50 mg(100/PK)	100 mg(100/PK)	200 mg(50/PK)	500 mg(50/PK)	500 mg(30/PK)	1000 mg(30/PK)
NH ₂	PP		730031		730033		
	Glass		730031G		730033G		
DMA	PP		730041		730043		
	Glass		730041G		730043G		
Diol	PP		730051		730053		
	Glass		730051G		730053G		
CN	PP		730061		730063		
	Glass		730061G		730063G		
NO ₂	PP				730143		730142
	Glass				730143G		730142G
Silica (SiOH)	PP		730071		730073	730070	730075
	Glass		730071G		730073G	730070G	730075G
PCA	PP		730480	730481	730482	730483	730484
	Glass		730480G	730481G	730482G	730483G	730484G
SA	PP		730076		730077		
	Glass		730076G		730077G		
PSA	PP		730460	730461	730462	730463	730464
	Glass		730460G	730461G	730462G	730463G	730464G
SB	PP		730078		730079		
	Glass		730078G		730079G		
Alox N	PP						730139
	Glass						730139G
Florisil	PP				730081		730082
	Glass				730081G		730082G
PA	PP					730007	730127
	Glass					730007G	730127G
HR-P	PP	730343		730108	730117	730111	730118
	Glass	730343G		730108G	730117G	730111G	730118G