

過濾頭 TITAN[®] Syringe Filter

4 mm 過濾頭

過濾樣品量 < 1ml

滯流容積 < 7.5µl

包裝數量 100/PK



材質	0.20µm	0.45µm
Nylon	42204-NN	44504-NN
PTFE	42204-NP	44504-NP
PVDF	42204-PV	44504-PV
Cellulose Acetate	42204-CA	44504-CA
Regenerated Cellulose	52204-RC	54504-RC

13 mm 過濾頭

過濾樣品量 < 10ml

滯流容積 < 25µl

包裝數量 100/PK



材質	0.20µm	0.45µm
Nylon	42213-NN	44513-NN
PTFE	42213-NP	44513-NP
PVDF	42213-PV	44513-PV
PP	42213-PP	44513-PP
Polysulfone	42213-PS	44513-PS
Cellulose Acetate	42213-CA	44513-CA
Regenerated Cellulose	52213-RC	54513-RC

25 mm 過濾頭

過濾樣品量 < 100ml

滯流容積 < 115µl

包裝數量 100/PK
*50/PK



材質	0.20µm	0.45µm
Nylon	42225-NN	44525-NN
PTFE	42225-NP	44525-NP
PVDF	42225-PV	44525-PV
PP (No Pre-Filter)	42225-PP	44525-PP
Polysulfone	42225-PS	44525-PS
Cellulose Acetate	42225-CA	44525-CA
材質	0.70µm	1.0µm
GMF	40725-GM	
PTFE		* 41025-NP
材質	1.2µm	3.1µm
GMF	41225-GM	42725-GM
Nylon	* 41225-NN	
材質	5.0µm	
Nylon	* 45025-NN	
PTFE	* 45025-NP	