



# 固相萃取裝置 Vacuum Manifold

## *Chromabond<sup>®</sup> Vacuum Manifold*



### *Vacuum Manifold and replacement parts*

目錄編號	規 格	包 裝
730150	Vacuum manifold complete consists of: glass cabinet with lid and lid gasket, removable needles on lower side of lid, vacuum gauge, control valve, valves and caps, variable rack, for 12 columns	1
730360	Vacuum manifold as described above, for 16 LV columns	1
730151	Vacuum manifold as described above, for 24 columns	1
730173	Glass cabinet without accessories, for 12 columns	1
730174	Glass cabinet without accessories, for 16 LV or 24 columns	1
730233	PP tank for vacuum manifold for 12 columns	2
730175	Lip with gasket, for 12 columns	1
730365	Lip with gasket, for 16 LV columns	1
730176	Lip with gasket, for 24 columns	1
730177	Gaskets for lid, for 12 columns	2
730178	Gaskets for lid, for 24 columns	2
730183	Luer fitting for lid, female	1
730184	Luer fitting for lid, male	1
730187	Drying attachment, for 12 columns	1
730188	Drying attachment, for 24 columns	1
730366	Collecting rack for 16 LV columns	1
730179	Vacuum gauge, complete with accessories	1
730152	Stainless steel needles	12
730154	Polypropylene needles	12
730185	Valve, plastic	12
730189.1	Valve, brass, tarnished	1
730189.12	Valve, brass, tarnished	12
730106	Stainless steel connectors	12
730243	Tubing adaptor for application of large sample volumes to polypropylene columns	4