

自動注射器樣品瓶 AutoSampler Vial

螺旋蓋樣品瓶 *Screw Cap Vial(Standard mouth 8-425 Thread) 1.5ml*



適用於各型式之自動注射系統，廠牌如：
HP, Varian, Shimadzu, Hitachi, Jasco, Perkin-Elmer, Waters, Beckman, LDC, Spectra-Physics, Thermal Separation 等。

目錄編號	規 格	包裝數量
200 250	透明樣品瓶 Screw Cap Vial, Clear Glass, 12 x 32mm, 1.5ml	100/pk
200 252	褐色樣品瓶 Screw Cap Vial, Amber Glass, 12 x 32mm, 1.5ml	100/pk
200 292E	有孔蓋 Open Top Cap, 8-425, Black PP	100/pk
200 292	有孔蓋 Open Top Cap, 8-425, Black PP	1000/pk
200 332	墊片 Septum Red Teflon/Silicone, 8mm, 0.060"	1000/pk
200 336	墊片 Septum Red Teflon/Silicone, 8mm, 0.060"	100/pk
200 368	墊片 Septum Teflon, 8mm, 0.010"	1000/pk
200 228E	平底內管 Insert 200µl, Flat Bottom, Sun51 Glass	100/pk
200 228	平底內管 Insert 200µl, Flat Bottom, Sun51 Glass	1000/pk
200 238	尖底內管 Insert 100µl, Conical, Sun51 Glass	100/pk
500 918	尖底內管附 PP 腳架 Insert 100µl, Mandrel Point, with PP Spring	100/pk

壓蓋式樣品瓶 *Crimp Top Vial 1.5ml*



適用於各型式之自動注射系統，廠牌如：
HP, Varian, Shimadzu, Hitachi, Jasco, Perkin-Elmer, Waters, Beckman, LDC, Spectra-Physics, Thermal Separation 等。

目錄編號	規 格	包裝數量
200 000	透明樣品瓶 Crimp Top Vial, Clear Glass, 12 x 32mm, 1.5ml	100/pk
200 002	褐色樣品瓶 Crimp Top Vial, Amber Glass, 12 x 32mm, 1.5ml	100/pk
200 104	鋁蓋 Aluminum Seal 11mm, Teflon /Rubber Septum	100/pk
200 100	鋁蓋 Aluminum Seal 11mm, Teflon /Rubber Septum	1000/pk
200 228E	平底內管 Insert 200µl, Flat Bottom, Sun51 Glass	100/pk
200 228	平底內管 Insert 200µl, Flat Bottom, Sun51 Glass	1000/pk
200 238	尖底內管 Insert 200µl, Conical, Sun51 Glass	100/pk
500 918	尖底內管附 PP 腳架 Insert 100µl, Mandrel Point, with PP Spring	100/pk

自動注射器樣品瓶 AutoSampler Vial

寬口螺旋蓋樣品瓶 *SUN1™ Screw Cap Vial (Wide mouth) 1.5ml*



HP 及 Varian Autosampler 專用

目錄編號	規 格	包裝數量
500 778	寬口透明樣品瓶 Screw Cap Vial, Clear Glass, 12 x 32mm, 1.5ml	100/pk
500 779	寬口褐色樣品瓶 Screw Cap Vial, Amber Glass, 12 x 32mm, 1.5ml	100/pk
500 790	黃色有孔蓋附墊片 Yellow Open Top Cap, with TFE/Sil Septum	100/pk
500 792	黃色有孔蓋附墊片 Yellow Open Top Cap, w/ TFE/Sil/TFE Septum	100/pk
501 790	綠色有孔蓋附墊片 Green Open Top Cap, with TFE/Sil Septum	100/pk
502 790	紅色有孔蓋附墊片 Red Open Top Cap, with TFE/Sil Septum	100/pk
503 790	黑色有孔蓋附墊片 Black Open Top Cap, with TFE/Sil Septum	100/pk
504 790	藍色有孔蓋附墊片 Blue Open Top Cap, with TFE/Sil Septum	100/pk
500 784	墊片 Septum Teflon/Silicone Rubber, 9mm, 0.040"	100/pk
500 822	墊片 Septum Teflon/Silicone Rubber 9mm, 0.040"	1000/pk
500 786	墊片 Septum Teflon/Silicone Rubber/Teflon, 9mm, 0.040"	100/pk
500 824	墊片 Septum Teflon/Silicone Rubber/Teflon, 9mm, 0.040"	1000/pk
200 668E	平底內管 Insert 350µl, Flat Bottom, Sun51 Glass	100/pk
200 668	平底內管 Insert 350µl, Flat Bottom, Sun51 Glass	500/pk
200 670	平底內管 Insert 350µl, Flat Bottom, Sun51 Glass	1000/pk
500 304	尖底內管 Insert 250µl, Conical, Sun33 Glass	100/pk
500 900	尖底內管附 PP 腳架 Insert 200µl, Mandrel Point, with PP Spring	100/pk

寬口螺旋蓋樣品瓶 *Screw Cap Vial (Wide mouth 10-425 Thread) 1.5ml*



Waters Alliance™ /2690

Autosampler 專用

目錄編號	規 格	包裝數量
500 306	寬口透明樣品瓶 Screw Cap Vial, Clear Glass, 12 x 32mm, 1.5ml	100/pk
500 322	寬口褐色樣品瓶 Screw Cap Vial, Amber Glass, 12 x 32mm, 1.5ml	100/pk
500 309	有孔蓋 Open Top Cap 10-425, Black PP	100/pk
500 894	墊片 Septum Tan Teflon/Silicone Rubber, 10mm, 0.060"	100/pk

※此樣品瓶所使用之內管與上列目錄編號 (200 668E, 200 668, 200 670, 500 304, 500 900) 相同。

宏濬儀器有限公司 <http://www.hplc.com.tw/>

見誠科技有限公司 <http://www.kkk.com.tw/>

自動注射器樣品瓶 AutoSampler Vial

塑膠螺旋蓋及壓蓋式樣品瓶 Polypropylene Vial 1.5ml



螺旋蓋塑膠樣品瓶



壓蓋式塑膠樣品瓶

螺旋蓋塑膠樣品瓶 Screw Cap PP Vial

目錄編號	規 格	包裝數量
200 268	塑膠樣品瓶 Screw Cap Vial, PP, 12 x 32mm, 250µl	100/pk
200 264	塑膠樣品瓶 Screw Cap Vial, PP, 12 x 32mm, 250µl	1000/pk
500 118	塑膠樣品瓶 Screw Cap Vial, PP, 12 x 32mm, 300µl	100/pk
500 114	塑膠樣品瓶 Screw Cap Vial, PP, 12 x 32mm, 300µl	1000/pk
500 124	塑膠樣品瓶 Screw Cap Vial, PP, 12 x 32mm, 500µl	100/pk
500 120	塑膠樣品瓶 Screw Cap Vial, PP, 12 x 32mm, 500µl	1000/pk
200 292E	有孔蓋 Open Top Cap, 8-425, Black PP	100/pk
200 292	有孔蓋 Open Top Cap, 8-425, Black PP	1000/pk
200 332	墊片 Septum Red Teflon/Silicone, 8mm, 0.060"	1000/pk
200 336	墊片 Septum Red Teflon/Silicone, 8mm, 0.060"	100/pk
200 368	墊片 Septum Teflon, 8mm, 0.010"	1000/pk

壓蓋式塑膠樣品瓶 Crimp Top PP Vial

目錄編號	規 格	包裝數量
200 050	塑膠樣品瓶 Crimp Top Vial, PP, 12 x 32mm, 250µl	100/pk
200 046	塑膠樣品瓶 Crimp Top Vial, PP, 12 x 32mm, 250µl	1000/pk
500 106	塑膠樣品瓶 Crimp Top Vial, PP, 12 x 32mm, 300µl	100/pk
500 102	塑膠樣品瓶 Crimp Top Vial, PP, 12 x 32mm, 300µl	1000/pk
500 112	塑膠樣品瓶 Crimp Top Vial, PP, 12 x 32mm, 500µl	100/pk
500 108	塑膠樣品瓶 Crimp Top Vial, PP, 12 x 32mm, 500µl	1000/pk
200 104	鋁蓋 Aluminum Seal 11mm, Teflon /Rubber Septum	100/pk
200 100	鋁蓋 Aluminum Seal 11mm, Teflon /Rubber Septum	1000/pk

自動注射器樣品瓶 AutoSampler Vial

螺旋蓋樣品瓶 *Screw Cap Vial (13-425 Thread) 4ml*



Waters 48-Position Carousel 專用

目錄編號	規 格	包裝數量
200 492	透明樣品瓶 Screw Cap Vial, Clear Glass, 15 x 45mm, 4ml	100/pk
200 494	褐色樣品瓶 Screw Cap Vial, Amber Glass, 15 x 45mm, 4ml	100/pk
500 576	有孔蓋 Open Top Cap, 13-425, Black PP	100/pk
200 594	墊片 Septum Red Teflon/Silicone, 12mm, 0.075"	100/pk
200 590	墊片 Septum Red Teflon/Silicone, 12mm, 0.075"	1000/pk
200 608	墊片 Septum Teflon, 12mm, 0.010"	1000/pk
200 756E	平底內管 Insert 750µl, Flat Bottom, Sun51 Glass	100/pk
200 756	平底內管 Insert 750µl, Flat Bottom, Sun51 Glass	1000/pk
200 744	塑膠尖底內管 Insert 250µl, Conical Bottom, Polypropylene	1000/pk
200 774	尖底內管 Insert 250µl, Conical Bottom, Sun51 Glass	100/pk
200 758	尖底內管 Insert 250µl, Mandrel Point, Glass with PP Spring	100/pk
200 778	金屬彈簧 Metal Spring for Inserts	100/pk



壓蓋式樣品瓶 *SUNVIALS™ 8 x 40 mm, 1ml*

Waters 96 & 48-Position Carousels 專用

目錄編號	規 格	包裝數量
500 602	透明樣品瓶(附蓋)SunVials, Clear Glass, 8 x 40mm, 1ml, w/PE cap	250/pk
500 634	透明樣品瓶(附蓋)SunVials, Clear Glass, 8 x 40mm, 1ml, w/PE cap	2000/pk
500 604	褐色樣品瓶(附蓋)SunVials, Amber Glass, 8 x 40mm, 1ml, w/PE Cap	250/pk
500 636	褐色樣品瓶(附蓋)SunVials, Amber Glass, 8 x 40mm, 1ml, w/PE Cap	2000/pk
500 765	塑膠樣品瓶(附蓋)ConSert PP SunVials, 8 x 40mm, 700µl, w/PE cap	250/pk

一般樣品瓶 Screw Cap Sampler Vial

螺旋蓋樣品瓶 8ml, 12ml, 16ml, 20ml, 24ml




目錄編號	規 格	包裝數量
200 330	無孔蓋 Solid Top Cap, 8-425, With TFE/ Rubber Liner	100/pk
200 568	無孔蓋 Solid Top Cap, 13-425, With TFE/ Rubber Liner	200/pk
203 498	透明樣品瓶組 Screw Cap Vial , Clear, 17 x 60mm, 8ml With Open Top Cap 15-425 & TFE/Silicone Septum	200/pk
203 499	褐色樣品瓶組 Screw Cap Vial , Clear, 17 x 60mm, 8ml With Open Top Cap 15-425 & TFE/Silicone Septum	200/pk
203 501	透明樣品瓶組 Screw Cap Vial , Clear, 19 x 65mm, 12ml With Open Top Cap 15-425 & TFE/Silicone Septum	144/pk
200 835	無孔蓋 Solid Top Cap, 15-425, With TFE/ Rubber Liner	200/pk
953 542	墊片 White PTFE/Silicone Septum, 15mm, 0.060"	100/pk
203 505	透明樣品瓶組 Screw Cap Vial , Clear, 21 x 70mm, 16ml With Open Top Cap 18-400 & TFE/Silicone Septum	144/pk
200 840	無孔蓋 Solid Top Cap, 18-400, With TFE/ Rubber Liner	200/pk
953 544	墊片 White PTFE/Silicone Septum, 18mm, 0.060"	100/pk
203 503	透明樣品瓶組 Screw Cap Vial , Clear, 28 x 57mm, 20ml With Open Top Cap 24-400 & TFE/Silicone Septum	144/pk
200 846	無孔蓋 Solid Top Cap, 24-400, With TFE/ Rubber Liner	200/pk
203 508	透明樣品瓶組 Screw Cap Vial , Clear, 23 x 85mm, 24ml With Open Top Cap 20-400 & TFE/Silicone Septum	200/pk
200 844	無孔蓋 Solid Top Cap, 20-400, With TFE/ Rubber Liner	100/pk
953 531	墊片 White PTFE/Silicone Septum, 20mm, 0.060"	100/pk

※此類樣品瓶可零買，但最少購買量為 50 個。

環境分析用樣品瓶 EPA Sampler Vial

螺旋蓋樣品瓶 40ml



目錄編號	規 格	包裝數量
200 900	透明樣品瓶 Screw Cap Vial , Clear, 27 x 95mm, 40ml	100/pk
200 902	褐色樣品瓶 Screw Cap Vial , Clear, 27 x 95mm, 40ml	100/pk
310 050	黑有孔蓋 Screw Open Top Cap, 24-414, Black, PP	100/pk
312 050	白有孔蓋 Screw Open Top Cap, 24-414, White, PP	100/pk
200 846	無孔蓋 Solid Top Cap, 24-400, With TFE/ Rubber Liner	200/pk
200 920	墊片 White PTFE/Silicone Septum, 22mm, 0.100"	100/pk
500 178	墊片 White PTFE/Silicone Septum, 22mm, 0.125"	100/pk
500 180	墊片 White PTFE/Silicone Septum, 22mm, 0.125"	1000/pk
200 906	透明樣品瓶組 Screw Cap Vial , Clear, 27 x 95mm, 40ml With Open Top Cap 24-414 & TFE/Silicone Septum 0.125"	100/pk
200 910	褐色樣品瓶組 Screw Cap Vial , Clear, 27 x 95mm, 40ml With Open Top Cap 24-414 & TFE/Silicone Septum 0.125"	100/pk

樣品瓶組 Assembled



目錄編號	規 格	包裝數量
200 912	透明樣品瓶組 Screw Cap Vial , Clear, 40ml, 27 x 95mm, with Black Cap/TFE/Sil Septum 0.125"	100/pk
200 913	透明樣品瓶組 Screw Cap Vial , Clear, 40ml, 27 x 95mm, with White Cap/TFE/Sil Septum 0.125"	100/pk
200 911	褐色樣品瓶組 Screw Cap Vial , Amber, 40ml, 27 x 95mm, with Black Cap/TFE/Sil Septum 0.125"	100/pk

微量反應樣品瓶 Micro-Reaction Vial



◎ 無刻度，附有孔蓋及 TFE/SIL 墊片

目錄編號	規 格	包裝數量
986281	微量反應瓶 React-Vial , 12 x 32mm, 0.1ml	12/pk
986282	微量反應瓶 React-Vial , 13 x 32mm, 0.3ml	12/pk
986284	微量反應瓶 React-Vial , 13 x 41mm, 1.0ml	12/pk
986288	微量反應瓶 React-Vial , 20 x 41mm, 2.0ml	12/pk
986287	微量反應瓶 React-Vial , 20 x 47mm, 3.0ml	12/pk
986289	微量反應瓶 React-Vial , 20 x 62mm, 5.0ml	12/pk

◎ 無刻度，附無孔蓋及 TFE/Rubber 墊片

目錄編號	規 格	包裝數量
986253	微量反應瓶 React-Vial , 13 x 32mm, 0.3ml	12/pk
986254	微量反應瓶 React-Vial , 13 x 41mm, 1.0ml	12/pk
986261	微量反應瓶 React-Vial , 20 x 41mm, 2.0ml	12/pk
986257	微量反應瓶 React-Vial , 20 x 47mm, 3.0ml	12/pk
986259	微量反應瓶 React-Vial , 20 x 62mm, 5.0ml	12/pk

◎ 有刻度，附有孔蓋及 TFE/SIL 墊片

目錄編號	規 格	包裝數量
986293	微量反應瓶 React-Vial , 13 x 32mm, 0.3ml	12/pk
986294	微量反應瓶 React-Vial , 13 x 41mm, 1.0ml	12/pk
986298	微量反應瓶 React-Vial , 20 x 41mm, 2.0ml	12/pk
986297	微量反應瓶 React-Vial , 20 x 47mm, 3.0ml	12/pk
986299	微量反應瓶 React-Vial , 20 x 62mm, 5.0ml	12/pk

◎ 有刻度，附無孔蓋及 TFE/Rubber 墊片

目錄編號	規 格	包裝數量
986273	微量反應瓶 React-Vial , 13 x 32mm, 0.3ml	12/pk
986274	微量反應瓶 React-Vial , 13 x 41mm, 1.0ml	12/pk
986278	微量反應瓶 React-Vial , 20 x 41mm, 2.0ml	12/pk
986277	微量反應瓶 React-Vial , 20 x 47mm, 3.0ml	12/pk
986279	微量反應瓶 React-Vial , 20 x 62mm, 5.0ml	12/pk

頂層空間取樣瓶 HeadSpace Vial



目錄編號	規 格	包裝數量
405 032	頂層空間取樣瓶 HeadSpaceVial , 22 x 38mm, 6ml, Rounded	100/pk
405 034	頂層空間取樣瓶 HeadSpaceVial , 22.5 x 45mm, 10ml, Rounded	100/pk
500 550	頂層空間取樣瓶 HeadSpaceVial , 22.5 x 75mm, 20ml, Rounded	100/pk
500 488	頂層空間取樣瓶 HeadSpaceVial , 22.5 x 45mm, 10ml, Square	100/pk
200 898	頂層空間取樣瓶 HeadSpaceVial , 22.5 x 75mm, 20ml, Square	100/pk
500 334	鋁蓋 Aluminum Cap 20mm, Standard Open-Top	1000/pk
200 998	鋁蓋 Aluminum Cap 20mm, Open-Top w/Retaining Ring	1000/pk
500 336	鋁蓋 Aluminum Cap 20mm, Pressure Release	1000/pk
405 038	鋁蓋 Aluminum Cap 20mm, Magnetic Seal*	1000/pk
400 910	塑膠蓋 Polyethylene Snap Plug 20mm	1000/pk
405 048	墊片 Teflon/Gray Butyl Septum, 20mm, 0.129"	100/pk
405 044	墊片 Teflon/Gray Butyl Septum, 20mm, 0.129"	1000/pk
500 772	墊片 Teflon/Red Rubber Septum, 20mm, 0.137"	1000/pk
200 952	墊片 Teflon/Silicone Septum, 20mm, 0.130"	100/pk
200 953	墊片 Teflon/Silicone Septum, 20mm, 0.130"	1000/pk
200 932	墊片塞 Gray Butyl Stopper, 20mm	1000/pk
405 061	墊片 Pharma-Fix Molded Butyl/Teflon Septum, 20mm, 0.130"	100/pk
405 063	墊片 Pharma-Fix Molded Butyl/Teflon Septum, 20mm, 0.130"	1000/pk
500 468	墊片塞 Butyl Rubber/PTFE, HeadLock Stopper, 20mm	1000/pk

* For use with Leap HS 500, CTC HS 500, Fisons HS 800.

Pharma-Fix™ Temperature range: -40°C to 120°C

儲存瓶 Graduated Bottle



玻璃儲存瓶



塑膠螺旋蓋(附管孔)



塑膠儲存瓶

目錄編號	規 格	包裝數量
201 682	玻璃儲存瓶 Clear Graduated Bottles 50ml, with 32mm PP Cap	1/ea
201 688	玻璃儲存瓶 Clear Graduated Bottles 100ml, with 45mm PP Cap	1/ea
201 696	玻璃儲存瓶 Clear Graduated Bottles 250ml, with 45mm PP Cap	1/ea
201 704	玻璃儲存瓶 Clear Graduated Bottles 500ml, with 45mm PP Cap	1/ea
201 720	玻璃儲存瓶 Clear Graduated Bottles 1000ml, with 45mm PP Cap	1/ea
201 720	玻璃儲存瓶 Clear Graduated Bottles 2000ml, with 45mm PP Cap	1/ea
998867	塑膠螺旋蓋 GL-45 Bottle Caps for 3/16"OD Tubing	1/ea
998755	塑膠螺旋蓋 GL-45 Bottle Caps for 1/8"OD Tubing	1/ea
998868	塑膠螺旋蓋 GL-45 Bottle Caps for 1/16"OD Tubing	1/ea
998800	塑膠螺旋蓋 GL-38 Bottle Caps for 1/8"OD Tubing	1/ea
998869	塑膠螺旋蓋 GL-38 Bottle Caps for 1/16"OD Tubing	1/ea
201 682P	塑膠儲存瓶 Clear Graduated Bottles 50ml, with 24mm PP Cap	24/cs
201 688P	塑膠儲存瓶 Clear Graduated Bottles 100ml, with 24mm PP Cap	24/cs
201 696P	塑膠儲存瓶 Clear Graduated Bottles 250ml, with 38mm PP Cap	12/cs
201 704P	塑膠儲存瓶 Clear Graduated Bottles 500ml, with 38mm PP Cap	12/cs
201 720P	塑膠儲存瓶 Clear Graduated Bottles 1000ml, with 55mm PP Cap	12/cs
201 720	塑膠儲存瓶 Clear Graduated Bottles 2000ml, with 55mm PP Cap	12/cs

玻璃樣品瓶 NEG Glass Vial



目錄編號	規 格	容 量	包 裝
透 明 clear			
H5-S1	15 x 45 mm	4 ml	100/pk
H5-S2	16.5 x 52 mm	5 ml	100/pk
H5-S3	19 x 55 mm	10 ml	100/pk
H5-S5	19 x 70 mm	12 ml	100/pk
H5-SV8	21 x 45 mm	8 ml	50/pk
H5-SV10	24 x 45 mm	10 ml	50/pk
H5-SV15	25 x 50 mm	15 ml	50/pk
H5-SV20	27 x 55 mm	20 ml	50/pk
H5-SV30	30 x 65 mm	30 ml	50/pk
H5-SV50	40 x 75 mm	50 ml	50/pk
H5-SV100	40 x 120 mm	100 ml	25/pk
H5-SV110	50 x 90 mm	110 ml	20/pk
H5-SV150	50 x 125 mm	150 ml	20/pk
褐 色 Amber			
H5-1S1	15 x 45 mm	4 ml	100/pk
H5-1S2	16.5 x 52 mm	5 ml	100/pk
H5-1S3	19 x 55 mm	10 ml	100/pk
H5-1S5	19 x 70 mm	12 ml	100/pk
H5-1SV8	21 x 45 mm	8 ml	50/pk
H5-1SV10	24 x 45 mm	10 ml	50/pk
H5-1SV15	25 x 50 mm	15 ml	50/pk
H5-1SV20	27 x 55 mm	20 ml	50/pk
H5-1SV30	30 x 65 mm	30 ml	50/pk
H5-1SV50	40 x 75 mm	50 ml	50/pk
H5-1SV100	40 x 120 mm	100 ml	25/pk
H5-1SV110	50 x 90 mm	110 ml	20/pk
H5-1SV150	50 x 125 mm	150 ml	20/pk

樣品瓶其他配件 Vial Other Products



目錄編號	規格	包裝數量
壓蓋器		
201 892	壓蓋器 Crimper 8mm	1/ea
200 246	壓蓋器 Crimper 11mm	1/ea
200 488	壓蓋器 Crimper 13mm	1/ea
200 972	壓蓋器 Crimper 20mm	1/ea
開蓋器		
500 440	開蓋器 Decrimper 8mm	1/ea
500 442	開蓋器 Decrimper 11mm	1/ea
200 248	開蓋器 Decapping Pliers 11mm	1/ea
500 444	開蓋器 Decrimper 13mm	1/ea
200 490	開蓋器 Decapping Pliers 13mm	1/ea
500 446	開蓋器 Decrimper 20mm	1/ea
200 974	開蓋器 Decapping Pliers 20mm	1/ea
樣品瓶架		
500 141	樣品瓶架 Vial Rack, 96 Position, White, PP, 8mm	1/ea
500 142	樣品瓶架 Vial Rack, 50 Position, White, PP, 12mm	1/ea
500 144	樣品瓶架 Vial Rack, 48 Position, White, PP 15mm	1/ea
200 442	樣品瓶架 Vial Rack, 100 Position, White, PP 12mm	1/ea
樣品瓶開關閥		
PS-614158	樣品瓶開關閥 Minimert Valves for 13-425 Thread Vials	12/ea
PS-614160	樣品瓶開關閥 Minimert Valves for 15-425 Thread Vials	12/ea
PS-614151	樣品瓶開關閥 Minimert Valves for 18-400 Thread Vials	12/ea
PS-614170	樣品瓶開關閥 Minimert Valves for 20-400 Thread Vials	12/ea
PS-614163	樣品瓶開關閥 Minimert Valves for 24-400 Thread Vials	12/ea
PS-614150	樣品瓶開關閥 Minimert Valves for Crimp Top 13mm ID Vials	12/ea
三角型攪拌子		
951698	三角型攪拌子 Magnetic Stirrers for 2.0, 3.0, 5.0ml Reaction Vials	6/pk
951699	三角型攪拌子 Magnetic Stirrers for 0.3, 1.0ml Reaction Vials	6/pk

塑膠針筒 Plastic Syringe



編號	容量	包裝
TOP-1	1 ml	100/pk
TOP-2.5	2.5 ml	100/pk
TOP-5	5 ml	100/pk
TOP-10	10 ml	100/pk
TOP-20	20 ml	50/pk
TOP-50	50 ml	25/pk

* 訂購時請告知是否需附針頭。