

液相層析儀 High Performance Liquid Chromatography

單流路等度液相層析儀 Isocratic analytical system



系統設備包含：

- 分析級輸送幫浦 – ECP2010 Analytical pump
- 紫外/可見光檢測器 – ECD2600 UV-VIS Detector
- 溶劑瓶架 – ECB2006 Box for bottles
- 樣品注入閥 – Analytical injection valve
- 層析控制軟體 – Clarity module for HPLC control
- 層析積分軟體 – Clarity Station for 1 instrument

四流路低壓梯度液相層析儀 Quaternary gradient analytical system



系統設備包含：

- 分析級輸送幫浦 – ECP2010 Analytical pump
- 紫外/可見光檢測器 – ECD2600 UV-VIS Detector
- 管柱加熱烘箱 – ECO2099 Column Oven
- 梯度控制系統及線上溶媒除氣裝置 – ECB2004 Gradient Box with Degasser
- 樣品注入閥 – Analytical injection valve
- 層析控制軟體 – Clarity module for HPLC control
- 層析積分軟體 – Clarity Station for 1 instrument

低壓梯度自動進樣液相層析儀 Gradient analytical system with autosampler



系統設備包含：

- 分析級輸送幫浦 – ECP2010 Analytical pump
- 紫外/可見光檢測器 – ECD2800 UV-VIS Detector
- 自動注射器 – L3320 Autosampler (for 2x48 vials)
- 管柱加熱烘箱 – ECO2080 Column Oven
- 梯度控制系統及線上溶媒除氣裝置 – ECB2004 Gradient Box with Degasser
- 樣品注入閥 – Analytical injection valve
- 層析控制軟體 – Clarity module for HPLC control
- 層析積分軟體 – Clarity Station for 1 instrument

液相層析儀 High Performance Liquid Chromatography

製備級等度液相層析儀 Isocratic Preparative system



系統設備包含：

- 製備級輸送幫浦 – ECP2100 Preparative Pump
- 紫外光檢測器 – Flash 14 DAD 400 Detector
- 溶劑瓶架 – ECB2005 Solvent box
- 樣品注入閥 – Analytical injection valve
- 層析控制軟體 – ECOMAC Software
- 層析積分軟體 – Peak-ABC Chromatography System

製備級梯度液相層析儀 Quaternary Gradient Preparative System



系統設備包含：

- 製備級輸送幫浦 – ECP2100 Preparative Pump
- 紫外/可見光檢測器 – Flash 14 DAD 600 Detector
- 梯度控制系統 – ECB2007 Gradient Box
- 樣品注入閥 – Analytical injection valve
- 層析控制軟體 – ECOMAC Software
- 層析積分軟體 – Peak-ABC Chromatography System

離子層析儀 Ion Chromatograph System



系統設備包含：

- 離子層析模組 – S150 Ion Chromatography Module
包含：電導度檢測器 Conductivity Detector
自動抑制器 Auto-Suppressor
管柱加熱器 Column Oven
- 非金屬輸送幫浦 – S1130 Isocratic Pump, PEEK
- 自動樣品進樣系統 – S 5300 Sample Injector System
- 離子層析管 – IC Column
- 層析積分軟體 – Data and Control System