

# 恆溫金屬浴器 Dry Bath



DB300-1



DB300-2



DB310-2

ChromTech® 熱蓋恆溫器規格：

產品型號	DB300-1	DB300-2	DB310-2
模塊數量	1 塊	2 塊	2 塊 (獨立分別控溫)
控溫範圍	室溫+5°C~150°C		
定時範圍	1min~99h59min (00:00 為連續)		
升溫速度	≤15 分鐘 (從 20°C 升至 150°C)		
溫度穩定性	≤±1°C (100~150°C) ≤±0.5°C (40~100°C)		
模塊溫度均勻性	≤±0.5°C		
顯示精度	0.1°C		
最大功率	200W	400W	
產品尺寸	W285×D225×H95mm		
重量	2.48 Kg	2.64 Kg	

DB300-1 / DB300-2 加熱塊規格表：

型號	規格	孔直徑	孔底形狀	模塊尺寸
DT01	6 mm x 42	6.5mm	圓底	95.5 x 76.5 x 50 mm
DT02	7 mm x 42	7.5mm	圓底	95.5 x 76.5 x 50 mm
DT03	10 mm x 20	10.5mm	圓底	95.5 x 76.5 x 50 mm
DT04	12 mm x 20	12.5mm	圓底	95.5 x 76.5 x 50 mm
DT05	13 mm x 20	13.5mm	圓底	95.5 x 76.5 x 50 mm
DT06	15 mm x 12 (7mL 離心管)	15.5mm	圓底	95.5 x 76.5 x 50 mm
DT07	16 mm x 12 (10mL/15mL 離心管)	16.5mm	圓底	95.5 x 76.5 x 50 mm
DT08	19 mm x 12	19.5mm	圓底	95.5 x 76.5 x 50 mm
DT09	20 mm x 6	20.5mm	圓底	95.5 x 76.5 x 50 mm
DT10	26 mm x 6	26.5mm	圓底	95.5 x 76.5 x 50 mm
DT11	28 mm x 4 (50mL 離心管)	28.5mm	平底	95.5 x 76.5 x 50 mm
DT12	40 mm x 2	40.5mm	圓底	95.5 x 76.5 x 50 mm
DT13	0.5mL x 42	8mm	錐底	95.5 x 76.5 x 50 mm
DT14	1.5mL x 24	10.8mm	錐底	95.5 x 76.5 x 50 mm
DT15	2.0mL x 24	10.8mm	圓底	95.5 x 76.5 x 50 mm
DT16	0.2mL x 48	6.1mm	錐底	95.5 x 76.5 x 50 mm

DB310-2 加熱塊規格表：

型號	規格	孔直徑	孔底形狀	模塊尺寸
DTD01	0.2mL x 48	6.1mm	錐底	95.5 x 69.3 x 50 mm
DTD02	0.5mL x 42	8mm	錐底	95.5 x 69.3 x 50 mm
DTD03	1.5mL x 24	10.8mm	錐底	95.5 x 69.3 x 50 mm
DTD04	2.0mL x 24	10.8mm	圓底	95.5 x 69.3 x 50 mm
DTD05	16 mm x 12 (10mL/15mL 離心管)	16.5mm	圓底	95.5 x 69.3 x 50 mm
DTD06	28 mm x 4 (50mL 離心管)	28.5mm	平底	95.5 x 69.3 x 50 mm