

往復式水浴震盪器 Water Bath Shaker

規格說明：



型號	WBS-10000
加熱容量	10L
溫度範圍	室溫(+5°C)~100°C
溫度精度	±0.5°C
定時範圍	1min~99h59m
震盪方式	往復式
振幅	26mm
震盪速度	20~200rpm
水槽尺寸	W250 x D304 x H200 mm
功率	1100W
電源	110V/60Hz or 220V/60Hz
外觀尺寸	W288 x D506 x H350 mm
淨重	10 Kg

恆溫水浴器 Thermostatic Water Bath

規格說明：



TWB-2000-1

TWB-2000-2

TWB-2000-3

TWB-5000

TWB-15000

型號	TWB-2000-1	TWB-2000-2	TWB-2000-3	TWB-5000	TWB-15000
容量	2L	2L x 2 孔	2L x 3 孔	5L	15L
溫度範圍	室溫加 5°C~100°C				
溫度精度	±0.5°C				
溫度顯示	LED 數字顯示 0.1°C				
溫度均溫性	±0.5°C				
定時範圍	1min~99h59min (0 為連續)				
功率	300W	600W	900W	700W	2200W
電源	110V/60Hz				220V/60Hz
外觀尺寸	W198 x D290 x H240 mm	W363 x D290 x H240 mm	W510 x D290 x H240 mm	W360 x D286 x H245 mm	W570 x D365 x H465 mm
淨重	2.3 Kg	5.1 Kg	6.7 Kg	5.1 Kg	13.6 Kg